

Container Labeling Format

last updated 7 August 2008

Make sure ALL vials and containers are properly labeled at ALL times. The label must have the following:

1. Contents
2. Concentration
3. Person responsible
4. Date

Blood Labeling Format:

1. Blood Type: A+, AB+, B+, O+, A-, AB-, B-, O-
2. Person responsible for the project (can use last names):
Smith
Jouben, etc.
3. Date Blood is drawn: MM / DD / YY

Example Label: **A+ - Leonard - 8 / 30 / 06**

This A+ blood was donated for use by Kaela on the 30th of August 2006.

It is very important that the donor name NEVER be connected with the sample in any records!